

# Datasheet

Cat#NBGN-100124

Version# RN5.9



---

<b>Product Name</b>	NebuChem™ Protein A/G Magnetic Beads (High Affinity)
<b>Size</b>	0.5ml; 5ml; Bulk
<b>Description</b>	NebuChem™ Protein A/G Magnetic Beads (High Affinity)(CAT#NBGN-100124) are ideal for small-scale purification of antibodies from plasma, serum, ascites, and cell culture.
<b>Storage and Stability</b>	This product is stable until the expiration date stated on the COA, when stored unopened at 2–8°C. Freezing is prohibited the product. Keep the magnetic beads in liquid suspension during storage and all handling steps. Drying will cause loss of binding capacity and result in reduced performance. Completely resuspend the beads before use. Observe sterile technique to avoid bacterial/fungal contamination.

---