

# Datasheet

Cat#NBGN-100117

Version# RN5.9



---

<b>Product Name</b>	NebuChem™ Anti-HisTag Affinity Resin
<b>Size</b>	1ml; 5ml; Bulk
<b>Description</b>	NebuChem™ Anti-HisTag Affinity Resin(CAT#NBGN-100117) is designed for purification of His-tagged protein. It is composed of agarose and conjugated THETM anti-His monoclonal antibody. His-tagged protein in the crude samples can bind to the anti-His monoclonal antibodies coupled with the agarose, and the nonspecific proteins can be eliminated by washing. As a result, the His-tagged protein will finally be eluted with high recovery. Moreover, this product can be used for immunoprecipitation of His-tagged protein.
<b>Storage</b>	Stored at 2-8 °C, Freezing is prohibited

---