

Datasheet

Cat# NBGN-100108

Version# RN5.9



Product Name	NebuChem™ Glutathione Resin
Size	10ml; 20ml; Bulk
Description	NebuChem™ Glutathione Resin(CAT#NBGN-100108) is an affinity chromatography medium designed for easy, one-step purification of recombinant glutathione S-transferase (GST) fusion proteins and other glutathione binding proteins expressed in E.coli, insect cells and mammalian cells. The recombinant GST fusion proteins can be purified directly from pre-treated cell lysate using Glutathione Resin. It is the excellent choice for high performance purifications.
Matrix	4% cross-linked agarose
Average Particle Size	90µm
Ligand	Glutathione
Binding Capacity	35-45mg GST-fusion protein (26kDa) /ml settled resin
Storage Solution	1 × PBS containing 20% ethanol
Storage	2-8 °C; Freezing is prohibited
